



# Multiplex SNaPshot for Body Fluid Identification

## Multiplex PCR

### Reagents Needed:

#### 4 X Primer Mix

AmpliTaq Gold® DNA Polymerase (Applied Biosystems, Foster City, CA)

Gold ST\*R 10 X Buffer (Promega, Madison, WI)

### 4 X Primer Mix for Multiplex PCR:

Target ID	Sequence (5'→3')	Conc. (uM)	Amplicon size (bp)	
SE1	cg17610929	TTG TTG ATA TGT TTT GAA TTA TTA AG	2.4	174
		ATA ACT TCC CTT ATC AAC ACC AAC	2.4	
SE2	Cg26763284-138d	TGA TTT ATA ATT ATT AGG GAG GGA AAT AG	0.8	105
		CCT AAA ACA ACC RAT TCC CAA C	1.6	
BL1	cg06379435	TTT ATT GGG GTA TTT TTA TTG GTT AG	8.0	157
		AAA ATA CAA CTT ACT CCT AAA CAC C	8.0	
BL3	cg08792630	TGT TTT AAG AGG ATG ATA AGG AA	2.4	220
		CCA CCT CAA TCC AAA CTA ACT ACA	2.4	
VF1	cg09765089-231d	TTG GTA GTT TTT GGA TTT TGG AG	2.4	137
		AAA CRT AAA ACR ACC CRA AC	19.2	
VF2	cg26079753-7d	TTT TGT GAG TGT GAG AGA TTT TTA AGA	1.6	176
		AAA ACC TCC AAA ACA AAA CCT CTA	1.6	
SA1	cg09652652-2d	GGG GAT TYG TTT YGT TAG GT	16.0	153
		CCA TTT CCC CCT TCC TAA AA	4.0	
MB	cg18069290	AGG GGY GAA GAG TAG GAA T	8.0	160
	cg09696411	ATA ATA AAA CAA CCA ATA CAA C	8.0	

### PCR Mixture:

PCR Component	Vol. (ul)
dH <sub>2</sub> O	~11.3
10 X Gold ST*R Buffer	2
4 X Primer Mix	5
AmpliTaq Gold (5 U/μL)	0.7 (3.5 U)
Bisulfite converted DNA	1 (~4)*
<b>Total</b>	<b>20</b>

### Thermal Cycling:

95°C for 11 minutes, then:

94°C for 20 seconds

56°C for 60 seconds

72°C for 30 seconds

for 34 cycles, then:

72°C for 7 minutes

4°C soak

\* Please be aware that you should not use too much volume of bisulfite converted DNA because, in our experience, it may cause PCR failure.

## Post-PCR Reaction

### Enzyme Purification of the PCR Product

#### Reagents Needed:

PCR product	5 $\mu$ L
ExoSAP-IT® (USB, Cleveland, OH)	2 $\mu$ L

#### Thermal Cycling:

37°C for 45 minutes
80°C for 15 minutes

### Multiplex SNaPshot

#### Reagents Needed:

5 X SBE Primer Mix  
5 X Sequencing buffer\_BigDye Termination (Applied Biosystems, Foster City, CA)  
SNaPshot™ Kit (Applied Biosystems, Foster City, CA)

#### 5 X SBE Primer Mix:

	Target ID	Sequence (5'→3')	Conc. (uM)	Length (nt)
SE1	cg17610929	(T) <sub>5</sub> CCG AAA CCC TCC CCA C	2.0	21
		(T) <sub>6</sub> CCA AAA CCC TCC CCA C	2.0	22
SE2	Cg26763284-138d	(T) <sub>6</sub> CGC GTA ACG ACT ATA AAA CCC TC	0.15	29
		(T) <sub>8</sub> CAC ATA ACA ACT ATA AAA CCC TC	0.5	31
BL1	cg06379435	(T) <sub>17</sub> CCR ATA AAA CCT CAA ACR TAA AAC	20.0	41
BL3	cg08792630	(T) <sub>21</sub> CCR TAA TAA CTT CTA CCT ATA AAT AAA CCC	3.0	51
VF1	cg09765089-231d	(T) <sub>34</sub> TCC CCA AAT AAC AAA CRA CRA AAA TC	4.5	60
VF2	cg26079753-7d	(T) <sub>44</sub> CRA TCA ACT ACT ATA AAA ACA CC	4.5	67
SA1	cg09652652-2d	(T) <sub>48</sub> CCA CGA ATA AAT AAC CAC GAT AAA AC	7.5	74
MB1	cg18069290	(T) <sub>61</sub> GAA CCA ACA ACA AAC RCC AC	3.0	81
MB2	cg09696411	(T) <sub>69</sub> CCT ATC CCR AAA CAA ACR C	15.0	88

#### SBE Reaction Mixture:

Reaction Component	Vol. (ul)
dH <sub>2</sub> O	~ 4
5 X SBE Primer Mix	2
5 X Sequencing Buffer	2
SNaPshot Reaction Mix	1
Purified PCR Product	> 1
Total	10

#### Thermal Cycling:

96°C for 10 seconds
50°C for 5 seconds
60°C for 30 seconds
for 25 cycles

\* Please keep the SNaPshot mixture on ice before putting it into the thermal cycler. Leaving the mixture at ambient temperature may result in a higher background because the SNaPshot kit does not support Hot Start PCR.

## Post-Single Base Extension

### Enzyme (SAP or CIP) Treatment

#### Reagents Needed:

SBE reaction product	10 $\mu$ L
SAP-Recombinant (USB, Cleveland, OH)	1 $\mu$ L

#### Thermal Cycling:

37°C for 45 minutes
80°C for 15 minutes

### Capillary Electrophoresis

#### Materials and Reagents Needed:

Dry heating block, water bath or thermal cycler  
3130 capillaries, 33 cm x 50  $\mu$ m (Applied Biosystems, Foster City, CA)  
Performance Optimized Polymer (POP4, Applied Biosystems, Foster City, CA)  
Matrix Standard Set DS-02 (dR110, dR6G, dTAMRA™, dROX™, LIZ® Dyes)  
(Applied Biosystems, Foster City, CA)  
Run Module GS STR POP4 (1 mL) E5  
GeneScan™ 120 LIZ™ Size Standard  
Hi-Di™ Formamide (Applied Biosystems, Foster City, CA)

#### Creating Matrix:

According to the ABI PRISM®SNaPshot™ Multiplex Kit protocol

#### Reagents Needed:

GeneScan™ 120 LIZ™ Size Standard	0.2 $\mu$ L
Hi-Di™ Formamide	10 $\mu$ L
SNaPshot product	1~2 $\mu$ L

#### Thermal Cycling:

95°C for 5 minutes
4°C soak


#### 3130 Data Collection Software:

Verify that you have chosen GeneScan Run module E5 and the DS-02 GeneScan Matrix Set.

Run prepared samples under the following conditions: injection time of 3 sec, electrophoresis voltage of 15 kV, collection time of 8 min, EP voltage of 15 kV and heat plate temperature of 60°C.

Detect and calculate peak heights with an analytical threshold of 100 rfu.

## Electropherograms



Representative electropherograms of body fluid identification using multiplex methylation SNaPshot; (A) semen, (B) blood, (C) vaginal fluid, (D) saliva, and (E) menstrual blood. SE1, SE2, BL1, BL3, VF1, VF2, SA1, MB1 and MB2 represent cg17610929, cg26763284-138d, cg06379435, cg08792630, cg09765089-231d, cg26079753-7d, cg09652652-2d, cg18069290 and cg09696411, respectively.